KSB-3 Complete Medium[®] Kit

Description

The KSB-3 Complete Medium[®] Kit has been developed for culture of mesenchymal stem cells (MSCs).

Content of the Product

Product	Catalog no.	Size	Qty	Storage	Shelf Life
KSB-3 Complete Medium® Kit	K3901	Kit	1	-	-
KSB-3 Basal Medium®	B1001	500mL	1	2-8°C*	12 months
KSB-3 Supplements®	S2901	2mL	1	-20°C*	12 months

* protect from light.

Important information

- KSB-3 Supplements® overnight at 4°C.

- KSB-3 Supplements[®] is stable up to at least one month at $2\sim 8^{\circ}$ C. Avoid refreeze.

- KSB-3 Complete Medium[®] is stable up to at least one month at $2 \sim 8^{\circ}$ C.

Medium preparation

- Add KSB-3 Supplements[®] to KSB-3 Basal Medium[®] (KSB-3 Complete Medium[®]).

- KSB-3 Complete Medium[®] requires supplementation with serum (commonly 5~10% Fetal Bovine Serum (FBS)).

- Antibiotic Supplement for Media (optional).

Guidelines for culture of mesenchymal stem cells (MSCs)

- We recommend refeeding the cultures every 2-3 days and passaging every 5-7 days.

- Do not let passaged MSCs become completely confluent, as it can reduce multipotency of MSCs.

- We recommend a seeding density of 3-6 x 10³ viable cells/cm².

Thawing of Cryopreserved MSCs

- **1.** Rapidly thaw frozen vial of cells in a 37°C water bath.
- 2. Transfer the cells into a 50mL conical tube.

 For every 1mL of cell suspension, add 10mL of KSB-3 Complete Medium[®].

4. Centrifuge cells at 200-400 x g for 5 minutes.

5. Resuspend cells in KSB-3 Complete Medium[®]. Take an aliquot from the cell suspension for cell counting.

6. Seed cells into culture vessels.

7. Incubate at 37° C in a humidified atmosphere containing 5% CO₂.

8. After 24hrs, discard spent medium and feed cells with fresh medium.

Refeed cells by complete medium exchanged every
2-3 days until cell passaging is needed.

Kangstem Biotech Co., Ltd.

Passing MSCs

- 1. Aspirate off spent medium.
- 2. Gently rinse cells with PBS and aspirate off.
- **3.** Add the 0.05~0.25% Trypsin-EDTA. Incubate vessels at 37° C.
- **4.** All the cells have detached, quickly proceed to the next step.

5. Add KSB-3 Complete Medium[®] to the vessel. Collect the cell suspension in a sterile 50mL conical tube.

6. Centrifuge cells at 200-400 x g for 5 minutes.

7. Remove supernatant.

8. Resuspend cells in KSB-3 Complete Medium[®]. Take an aliquot from the cell suspension for cell counting.

9. Seed cells into culture vessels.

10. Incubate at 37° C in a humidified atmosphere containing 5% CO₂.

Refeed cells by complete medium exchanged every
2-3 days until cell passaging is needed.

Cryopreservation of cells

1. Prepare cryopreservation solution by supplementing the KSB-3 Complete Medium[®] with 10% dimethyl sulfoxide (DMSO).

2. Harvest cells as described in the passing protocol.

3. Pellet cells by centrifugation at 200-400 x g for 5 minutes, resuspend cells in the cryopreservation solution.

4. Take an aliquot from the cell suspension for cell counting.

5. Transfer to cryogenic vials. 1°C freezing container and place in a -80°C freezer overnight.

6. Transfer cryogenic vials to the liquid nitrogen for long-term storage.

For Research Use Only. Not for use in diagnostic procedure.